

Table S1. Identity of internal standards and reference standards used in LC-MS/MS analysis. Lipid standards for major phospholipid classes were used in MS analysis (Avanti Polar Lipids, Inc., Alabaster, AL): phosphatidic acid (PA), phosphatidylcholine (PC), phosphatidylethanolamine (PE), phosphatidylglycerol (PG), phosphatidylinositol (PI), phosphatidylserine (PS), and cardiolipin (CL).

Lipid Class	Internal Standard	Reference Standard
PA	PA(16:0 D31-18:1)	PA(18:1-18:1)
PC	PC(16:0 D31-18:1)	PC(18:1-18:1)
PE	PE(16:0 D31-18:1)	PE(18:1-18:1)
PG	PG(16:0 D31-18:1)	PG(18:1-18:1)
PI	PI(16:0 D31-18:1)	PI(18:1-18:1)
PS	PS(16:0 D31-18:1)	PS(18:1-18:1)
CL	CL(14:0-14:0-14:0-14:0)	CL(18:2-18:2-18:2-18:2)