

Supplementary Table 1. Primers for FCGR3A-Sanger, FCGR3B-Sanger and FCGR3(A+B)-Sanger assays

Primer *	Function	Direction	Assay(s)	Sequence (5'-3')
1	PCR	Forward	FCGR3A-Sanger, FCGR3B-Sanger	CCC TTC ACA AAG CTC TGC ACT
4	PCR	Reverse	FCGR3A-Sanger	ATT CTG GAG GCT GGT GCT ACA
4	PCR	Reverse	FCGR3B-Sanger	ATT CTG GAG GCT GGT GCT ATG
2	Sequencing	Forward	FCGR3A-Sanger, FCGR3B-Sanger, FCGR3(A+B)-Sanger	GTT CAA GGA GGA AGA CCC TAT TC
3	Sequencing	Reverse	FCGR3A-Sanger, FCGR3B-Sanger, FCGR3(A+B)-Sanger	TGC AAA CCT ACC CTG CAA TC
2	PCR	Forward	FCGR3(A+B)-Sanger	GTT CAA GGA GGA AGA CCC TAT TC
5	PCR	Reverse	FCGR3(A+B)-Sanger	GCG TAA GGA AAG ACA GAG GTA TT

* refers to Supplementary Figures 1 and 2