

TABLE E-1 Settings for MRI Technique

	Sagittal		Axial	
	T1-Weighted	T2-Weighted	T1-Weighted	T2-Weighted
Repetition time (<i>ms</i>)	501-777	3000-3030	576-877	3000-5000
Time to echo (<i>ms</i>)	9.2	86	8.7	83-85
Field of view (<i>mm</i>)	300 × 300	300 × 300	260 × 260	260 × 260
Slice thickness (<i>mm</i>)	4	4	4	4

TABLE E-2 Settings for CT Imaging Technique

Slice thickness	0.6-2.0 mm
Pitch	0.6
Rotation time	1.5 sec
Peak kilovoltage	140 kV
Milliamperere-second	450 mAs
Bone filter	“D”
Axial reconstructions	2 mm thick
Coronal reconstructions	2 mm thick
Sagittal reconstructions	2 mm thick