

Supplemental Digital Content 2**CT Protocol chest****Oral Contrast:**

N

Pre Contrast:

O

Post Contrast:

N

Delay:

O

Bolus Tracking / Smart Prep:

No

Kernel / Algorithm:

B40 / B60

Scan Range:

Apices to Kidney

Positioning:**Supine/ * Pro ne****Recon #1:**

1mm x 1mm, B40, Mediastinum

Recon #2:

1mm x 1mm, B60, Lung

Recon #3:

Axial, MIP, 10mm x 5mm, B40, Lung

Recon #4:

Coronal, MPR, 2mm x 2mm, B40, Lung

Recon #5:

Sagittal, MPR, 2mm x 2mm, B40, Lung

Recon #6:

Expiration 1mm x 1mm B60, Lung (Optional)

Recon #7:

* Pro ne 1mm x 1mm, B60, Lung (Optional)

Scan Range: Carina to Diaphragm (Optional)**Technique:****kV:**

120

mA:**DFOV:**

Fit to Patient

Window:**Mediastinum:****WW=** 450**WL=** 45**Lung:****WW=** 1500**WL=** -700**Comments:****RIS Codes:**

CTC

CT Chest

