

Supplemental Digital Content 3: Sensitivity analysis - percentile ranks for systolic pressure variation & pulse pressure variation value distributions, cases with anesthesia induction end 09:00-13:00

	Percentile										
	1 st	2.5 th	5 th	10 th	25 th	50 th	75 th	90 th	95 th	97.5 th	99 th
SPV, Manually-Entered (mmHg) N = 739	2.5	3.0	3.5	4.0	5.0	6.0	8.0	10.0	11.0	13.0	14.0
SPV, Automated Recorded (mmHg) N = 214	2.0	2.1	2.3	2.9	3.8	4.7	6.1	7.9	9.9	13.1	15.7
PPV, Automated Recorded (%) N = 214	2.0	2.0	2.5	3.0	5.0	7.0	10.0	12.5	18.0	21.0	24.0

SPV = systolic pressure variation, PPV = pulse pressure variation